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- Balows, A., W. J. Hausler, Jr., K. L. Hemann, H. D. Isenberg, and J. H. Shadomy. 1991. Manual of Clinical Microbiology. 5th Edition. American Society for Microbiology, Washington D. C.
- Blazevic, D. J. and G. M. Ederer. 1975. Principals of Biochemical Tests in Diagnostic Microbiology. John Wiley & Sons, New York, N. Y.
- Cemey, G. 1978. Studies on the Aminopeptidase Test for the Distinction of Gram-negative from Gram-positive Bacteria. Europ. J. Appl. Microbiol. 5:115-122.
- Cowan, S. T. 1974. Cowan and Steel's Manual for the Identification of Medical Bacteria (2nd ed). Cambridge University Press, London.
- Eriquez, L. A., A. P. Jones, and N. E. Hodinka. 1992. Enzymatic Test System for the Rapid Identification of Enterobacteriaceae. Abstract C5. 92nd Annual Meeting, American Society for Microbiology, Washington DC.
- Eriquez, L. A., N. E. Hodinka, and K. A. Hornsby. 1993. A Survey of Clinically Relevant Enterobacteriaceae and Other Gram-negative Rods for β -Glucuronidase Activity. Abstract C294. 93rd Annual Meeting. American Society for Microbiology, Washington D. C.
- Ewing, W. E. 1986. Edwards and Ewing's Identification of Enterobacteriaceae. 4th Edition. Elsevier Science Publishing Company, Inc., New York, N. Y.
- Farmer III, J. J., B. R. Davis, F. W. Hickman-Brenner, A. McWhorter, G. P. Huntley-Carter, M. A. Asbury, C. Riddle, H. G. Wathen-Grady, C. Elias, G. R. Fanning, A. G. Steigerwalt, C. M. O'Hara, G. K. Morris, P. B. Smith, and D. J. Brenner. 1985. Biochemical Identification of New Species and Biogroups of Enterobacteriaceae Isolated from Clinical Specimens. J. Clinical Microbiol. 21:46-76.
- Farmer III, J. J., and M. T. Kelly. 1991. Enterobacteriaceae. In: Manual of Clinical Microbiology, 5th Edition., A. Balows Ed. in Chief. pp. 360-383. American Society for Microbiology, Washington, D. C.
- Finegold, S. M., W. J. Martin, and E. G. Scott. 1978. Bailey and Scott's Diagnostic Microbiology (5th ed.). C. V. Mosby Co., St. Louis, MO.
- Gilardi, G.L., S. Hirsch and M. Mandel. 1975. Characteristics of Yellow-pigmented Nonfermentative Bacilli (groups Ve-1 and Ve-2) Encountered in Clinical Bacteriology. J. Clin. Microbiol. 1:384-389.

- Gilardi, G. L. 1978. Glucose Nonfermenting Gram-negative Bacteria in Clinical Microbiology, CRC Press, West Palm Beach, FL.
- Guilbert, G. G. 1970. Determination of Enzymes. In: Methods of Enzymatic Analysis. p. 43-51. Pergamon Press, New York, NY.
- Holt, J. G. III and N. R. Kreig IV. 1984. Bergey's Manual of Determinative Bacteriology. Vol. 1. Williams and Wilkens, Baltimore, MD.
- Humble, W. M. A. King, and I. Phillips. 1977. Api Zym. A Single Rapid System for the Detection of Bacterial Enzymes. J. Clin. Pathology. 30:275-277.
- Isenberg, H. D. 1992. Clinical Microbiology Procedures Handbook. Vol. 1. American Society for Microbiology, Washington, DC.
- Jerris, R., A. Kabant, S. Peroulas, T. Lee, K. A. Homsby, and D. Lockhart. 1993. Comparison of the IDS RapID E Plus and the API 20E for Identification of Routine Clinical Isolates. Abstract C304. 93rd Annual Meeting, American Society for Microbiology, Washington, DC.
- Kilian, M. and P. Bulow. 1976. Rapid Diagnosis of Enterobacteriaceae I. Detection of Bacterial Glycosidases. Acta Pathol. Microbiol. Scand. Sect. B. 84:245-251.
- Kitch, T., M. R. Jacobs, and P. C. Applebaum. 1993. Evaluation of the 4-hour E Plus System for the Identification of 104 Enterobacteriaceae and Gram-negative Non-fermentors. Abstract C303. 93rd Annual Meeting. American Society for Microbiology, Washington DC.
- MacFaddin, J. F. 1980. Biochemical Tests for Identification of Medical Bacteria (2nd ed.). Williams and Wilkens, Baltimore, MD.
- Mulczyk, M. and A. Szewczuk. 1970. Pyrrolidonyl Peptidase in Bacteria. A New Colorimetric Test for Differentiation of Enterobacteriaceae. J. Gen. Microbiol. 61:9-13.
- Nagatsu, T., M. Hino, H. Fuyamada, T. Hayakowa, S. Sakakibara, Y. Nakagawa, and T. Takemoto. 1976. New Chromogenic Substrates for X-prolyl Dipeptidyl- aminopeptidase, Anal. Biochem. 74:466-476.
- Nord, C. E., A. A. Lindberg, and A. Dahlback. 1975. Four Hour Tests for the Identification of Enterobacteriaceae. Med. Microbiol. Immunol. 161:231-238.
- Peterson, E. H. and E. J. Hsu. 1978. Rapid Detection of Gram-negative Bacteria by Aminopeptidase Profiles. J. Food. Sci. 43:1853-1856.
- Schreckenberger, P., M. Montero, and N. Heldt. 1993. Evaluation of the Rapid E Plus Panel for Identification of Enterobacteriaceae. Abstract C309. 93rd Annual Meeting. American Society for Microbiology, Washington, DC.
- Waston, R. R. 1976. Substrate Specificities of Aminopeptidases. A Specific Method for Microbial Identification. Methods in Microbiol. 9:1-14.
- Weaver, R. E., D. G. Holiis, W. A. Clark and P. Riley. 1983. The Identification of Unusual Pathogenic Gram-negative Bacteria. Centers for Disease Control, Atlanta, GA.
- Westley, J. R., P. J. Anderson, V. A. Close, B. Halpem, and E. M. Lederberg. 1967. Aminopeptidase Profiles of Various Bacteria. Appl. Microbiol. 15:822-825.

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