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1. Arai, T. M., M. Otake, S. Enomoto, S. Goto, and S. Kuwahara, 1970. Determination of *Pseudomonas aeruginosa* by biochemical test methods. II. Acylamidase test, a modified biochemical test for the identification of *Pseudomonas aeruginosa*. Jpn. J. Microbiol. 14:279-284.
2. Baumann, P., M. Doudoroff, and R. Y. Stanier, 1968. A study of the *Moraxella* group II. Oxidase-negative species (genus *Acinetobacter*). J. Bacteriol. 95:1520-1541.
3. Blazevic, D. J. and G. M. Ederer, 1975. Principles of biochemical tests in diagnostic microbiology. John Wiley & Sons, New York, NY.
4. Cerney, G., 1978. Studies on the aminopeptidase test for the distinction of Gram-negative from Gram-positive bacteria. Europ. J. Appl. Microbiol. 5:113-122.
5. Clark, W.A., D. G. Hollis, R. E. Weaver, and P. Riley, 1984. The identification of unusual pathogenic Gram-negative bacteria, Centers for Disease Control. Atlanta, GA.
6. Cowan, S. T. 1974. Cowan and Steel's manual for the identification of medical bacteria, 2nd Edition, Cambridge

University Press, London.

7. Eriqez, L. A., A. P. Jones, and N. E. Hodinka. 1991. Enzymatic test system for the rapid identification of non-fermentative and selected glucose-fermenting Gram-negative rods. Abstract C217. Annual Meeting, American Society for Microbiology, Washington, DC.
8. Finegold, S. M., W. J. Martin, and E. G. Scott. 1978. Bailey and Scott's Diagnostic Microbiology, 5th Ed. C. V. Mosby Co., St. Louis Mo.
9. Gilardi, G. L. 1972. Infrequently encountered *Pseudomonas* species causing infections in humans. Ann. Intern. Med. 72:211-215.
10. Gilardi, G. L. 1976. *Pseudomonas* species in clinical microbiology. Mt. Sinai J. Med. 43:710-726
11. Gilardi, G. L., S. Hirschi, and M. Mandel. 1975. Characteristics of yellow-pigmented nonfermentative bacilli (groups Ve-1 and Ve-2) encountered in clinical bacteriology. J. Clin. Microbiol. 1:384-389.
12. Gilardi, G. L. 1978. Glucose nonfermenting Gram-negative bacteria in clinical microbiology. CRC Press, West Palm Beach, FL.
13. Guilbert, G. G. 1970. Determination of Enzymes. In: Methods of Enzymatic Analysis, p. 43-51. Pergamon Press, New York, NY.
14. Henriksen, S. D. 1973. *Moraxella*, *Acinetobacter*, and the *Mimeia*, Bacteriol. Rev. 37:522-561.
15. Holt, J., N. R. Krieg, P. H. A. Sneath, J. T. Stanley, and S. T. Williams. 1994. Bergey's Manual of Determinative Bacteriology, 9th ed. Williams and Wilkins, Baltimore, MD.
16. Hugh, R. 1970. A practical approach to the identification of certain nonfermentative Gram-negative rods encountered in clinical specimens. J. Conf. Public Health Lab. Directors. 28:168-187.
17. Humble, W. M., A. King, and I. Phillips. 1977. Api Zym. A single rapid system for the detection of bacterial enzymes. J. Clin. Pathol. 30:275-277.
18. Kilian, M. and P. Bulow. 1976. Rapid diagnosis of Enterobacteriaceae. I. Detection of bacterial glycosidases. Acta Pathol. Microbiol. Scan. Sect. B. 84:245-251.
19. Kiska, D. L., A. Kerr, M. C. Jones, J. A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P. H. Gilligan. 1995. Evaluation of four commercial identification systems for glucose non-fermenters from cystic fibrosis patients. Abstract C312. Annual Meeting, American Society for Microbiology, Washington, DC.
20. Kiska, D. L., A. Kerr, M. C. Jones, J. A. Caracciolo, B. Eskridge, M. Jordan, S. Miller, D. Hughes, N. King, and P. H. Gilligan. 1996. Accuracy of four commercial systems for the identification of *Burkholderia cepacia* and other Gram-negative nonfermenting bacilli recovered from patients with cystic fibrosis. J. Clinical Microbiol. 34:886-891.
21. Kitch, T., M. R. Jacobs, and P. C. Appelbaum. 1991. Evaluation of the 4hr NF Plus method for identification of 86 Gram-negative non-fermentative rods. Abstract C215. Annual Meeting, American Society for Microbiology, Washington, DC.
22. Kitch, T., M. R. Jacobs, and P. C. Appelbaum. 1992. Evaluation of the 4-hour RapID NF Plus method for identification of 345 Gram-negative non-fermentative rods. J. Clin. Microbiol. 30:1267-1270.
23. MacFadden, J. F. 1980. Biochemical tests for identification of medical bacteria, 2nd Ed. Williams and Wilkins, Baltimore, MD.
24. Mulczyk, M. and A. Szewczuk. 1970. Pyrrolidonyl peptidase in bacteria. A new colorimetric test for differentiation of Enterobacteriaceae. J. Gen. Microbiol. 61:9-13.
25. Murray, P.R., E. J. Baron, M. A. Pfaller, F. C. Tenover, and R. H. Tenover, Eds. 1995. Manual of Clinical Microbiology, 6th Edition, American Society for Microbiology, Washington, DC.
26. Nagatsu, T., M. Hino, H. Fuyamada, T. Hayakawa, S. Sakakibara, Y. Nakagawa, and T. Takemoto. 1976. New chromogenic substrates for x-prolyl dipeptidyl-aminopeptidase. Anal. Biochem. 74:466-476.
27. Nord, C. E., A.A. Linberg, and A. Dahlback. 1975. Four hour tests for the identification of Enterobacteriaceae. Med. Microbiol. Immunol. 161:231-238.
28. Peterson, E. H. and E. J. Hsu. 1978. Rapid detection of

- Gram-negative bacteria by aminopeptidase profiles. *J. Food Sci.* 43:1853:1856.
29. Pickett, J. M., D. G. Hollis, and E. J. Bottone. 1991. Miscellaneous Gram-negative bacilli, p. 410-428. In: *Manual of Clinical Microbiology*, 5th Ed. American Society for Microbiology, Washington, DC.
30. Stanier, R.Y., N. J. Palleroni, and M. Doudoroff. 1966. The aerobic pseudomonads: a taxonomic study. *J. Gen. Microbiol.* 43:159-271.
31. Watson, R. R. 1976. Substrate specificities of aminopeptidases. A specific method for microbial identification. *Methods In Microbiol.* 9:1-14.
32. Westley, J. R., P. J. Anderson, V. A. Close, B. Halpern, and E. M. Lederberg. 1967. Aminopeptidase profiles of various bacteria. *Appl. Microbiol.* 15:822-825

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