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1. Baron, E. J., L. R. Peterson and S. M. Finegold, Eds. 1994. Bailey & Scott's Diagnostic Microbiology, 9th ed. Mosby, St. Louis.
2. Cartwright, C. P., F. Stock, P. M. Kruczak-Filipov and V. J. Gill. 1993. Rapid method for presumptive identification of *Corynebacterium jeikeium*. J. Clin. Microbiol. 31:3320-3322.
3. Coyle, M. B. and B. A. Lipsky. 1990. Coryneform bacteria in infectious disease: clinical and laboratory aspects. Clin. Microbiol. Rev. 3:227-246.
4. Eriquez, L. A. and J. K. Marler. 1996. Rapid differentiation and identification of corynebacteria and other selected gram-positive bacilli using single-substrate enzymatic tests. Abstract C-393, Abstracts of the 96th General Meeting of the American Society for Microbiology, Washington, DC.
5. Freney, J., et. al. 1991. Evaluation of API Coryne in comparison with conventional methods for identifying Coryneform bacteria. J. Clin. Microbiol, 29:38-41.
7. Funke, G., A. Von Graevenitz, J. E. Clarridge III, and K. A. Bernard. 1997. Clinical Microbiology of Coryneform Bacteria, Clin. Microbiol. Rev. 10:125-159.
8. Gavin, S. E., et. al. 1992. Evaluation of the Rapid CORYNE identification system for *Corynebacterium* species and other coryneforms. J. Clin. Microbiol. 30:1692-1695.
9. George, M. J. 1995. Clinical significance and characterization of *Corynebacterium* species. Clin. Microbiol. Newsletter. 17:177-179.
10. Grasmick, A. E. and D. A. Bruckner. 1987. Comparison of rapid identification methods and conventional substrates for identification of *Corynebacterium* Group JK isolates. J. Clin. Microbiol. 25:1111-1112.
11. Guilbert, G. G. 1970. Determination of Enzymes. In: Methods of Enzymatic Analysis, pp. 43-51. Pergamon Press, New York, NY.
12. Hollis, D. G., F. O. Sottnek, W. J. Brown and R. E. Weaver. 1980. Use of the rapid fermentation test in determining carbohydrate reactions of fastidious bacteria in clinical laboratories. J. Clin. Microbiol. 12:620-623.
13. Hudspeth, M. et. al. 1997. Premarket evaluation prototype panel for identification of *Corynebacterium* species and Gram-positive rods. Abstract C459, Abstracts of the 97th General Meeting, American Society for Microbiology, Washington DC.
14. Isenberg, H. D. 1992. Clinical Microbiology Procedures Handbook. Vol. 1. American Society for Microbiology, Washington, DC.
16. MacFaddin, J. F. 1980. Biochemical Tests for identification of Medical Bacteria. 2nd Ed. Williams and Wilkins, Baltimore, MD.
17. Murray, P. E. J. Baron, M. A. Pfaller, F. C. Tenover and R. H. Tenover, Eds. 1995. Manual of Clinical Microbiology, 6th ed. American Society for Microbiology, Washington DC.
18. O'Leary, W. M. 1989. Practical Handbook of Microbiology. CRC Press, Boca Raton, FL.
19. Sneath, P. H., N. S. Mair, M. E. Sharpe and J. G. Holt, Eds. 1986. Bergy's Manual of Systematic Bacteriology. Williams and Wilkins, Baltimore, MD.
20. Watson, R. R. 1976. Substrate specificities of aminopeptidases: a specific method for microbial differentiation, pp. 1-14. In: Methods in microbiology, vol. 9. J. R. Norris and D. W. Ribbons, Eds. Academic Press, New York, NY.
21. Watson, R. R. 1976. Substrate specificities of aminopeptidases. A specific method for microbial identification. Methods in Microbiol. 9:1-14

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